

**Aurum Ceramic / Classic Scanning Facilities Across Canada!**

# PROCERA® AllCeram

**Superb All-Ceramic Esthetics. . . Prep and Cement Like A PFM**



**Natural Beauty.** Vital translucent esthetics and perfect shade blending without opacity. Dark or discoloured margins are eliminated. Ideal for covering gold or metal posts.

**Exceptional Strength.** Densely sintered zirconia oxide coping offers a biaxial strength (1121 MPa) far exceeding its milled competition (930 MPa) or Alumina (687 MPa). High initial strength helps prevent cracks, fatigue, stress and corrosion.

**Precision Fit.** Coping milled directly on the die using CAD/CAM technology.

**Traditional PFM Prep and Cementation.** Only basic chamfer preparation required. Combines esthetic freedom of using all-ceramic single crowns with the convenience of using your favourite PFM cement.

**Guaranteed.** Over 1,000,000 patients have been treated world-wide with Procera AllCeram! Each Procera AllCeram restoration comes with our 5-Year Warranty against breakage or defects.

## **Available Coast-to-Coast from Aurum Ceramic/Classic**

Procera AllCeram scanners have been placed in our Victoria, Vancouver, Calgary, Toronto and Ottawa laboratories. This allows us to deliver the quick turnarounds you expect while maintaining the superb craftsmanship you demand.


Aurum Ceramic/Classic has the knowledge, experience, proven technical abilities and creativity to deliver unsurpassed Procera AllCeram results on every case – every time. For more information, call your closest Aurum Ceramic/Classic laboratory TOLL FREE.




*Restorations fabricated by Aurum Ceramic/Classic.*

## **Procera® AllCeram Preparation Guide**

### **Posterior Preparation**




Finish line should be chamfer 0.5 to 0.7 mm subgingivally. Maintain relatively level topography when preparing occlusal surface.




Eliminate undercuts and sharp edges. Avoid steep slopes and sharp grooves.

### **Anterior Preparation**



Shape lingual surface with diamond bur to provide space for crown contour and retention.



Fine finish lines enable scanner to register all aspects of die models for optimal marginal fit.

