

# ECLIPSE™

*Taking a processed baseplate and building on it!*

## Monomer-free!

- Does not contain methyl, ethyl, propyl or butyl methacrylate monomers.

## A revolutionary new concept in denture fabrication

- No flasking, investing or boiling out. No processing errors or changes in vertical dimension.
- Light-cured denture base materials handle like wax.
- Final baseplate formed directly on Master Model. Used for occlusal records and try-in.
- Eclipse material and teeth easily adjusted chairside.
- Final processing accomplished by technician with high-energy light-curing system.



## Fit at Try-in equals fit at delivery

- Baseplate used for occlusal rims becomes part of final denture.
- Dentist and patient can evaluate function and aesthetics at the occlusal and try-in appointments.
- Enables patients to feel final fit prior to delivery, increasing final case acceptance.

## Clinically proven


- Multi-site clinical study conducted with over 100 patients and up to two years data.
- Superior fit and comfort. High strength.

## Easy to implement in your practice

- Customized aesthetics and easy characterization.
- Use acrylic or porcelain teeth.
- Fully repairable and relineable.

## Build Your Denture Practice With Aurum Ceramic

For more information on monomer-free Eclipse dentures, call TOLL FREE.



## Simplified indirect technique



Technician forms final baseplate directly on Master Model.



After Occlusal appointment, teeth set directly on baseplate.



Denture base and teeth easily adjusted chairside at Try-in appointment.