

5X The Choice

The #1 choice for over 10 years in Germany.*



- Futar®
- Futar® D
- Futar® D Fast
- Futar® D Slow
- Futar® Scan

PVS based bite registration materials.

* As reported in GfK

1-877-532-2123



In an independent clinical study conducted by a leading research publication, THE DENTAL ADVISOR, Futar® Bite Registration Materials received a 4 1/2 plus rating out of a possible 5 and was given the Editors' Choice Award. To read the full clinical evaluation article visit www.dentaladvisor.com.

Precision & Accuracy

When choosing a bite registration material,

would you prefer a high final hardness and comfortable working time? Or is speed the vital factor? Should your material be "scanable"? Whatever your requirements may be, you have a choice within the Futar® brand bite registration materials.

Simply choose the final hardness and working time you require, and benefit from the advantages of the Futar® brand's maximum precision and comfort. And you can take comfort in the fact that Futar® brand bite registration materials are international market leaders in Germany.



5X
The Advantages
Futar® Brand



- Extremely high final hardness for precise restorations
- Comfortable working times
- Can be conveniently milled and cut with a scalpel
- Highly thixotropic, will not run into the interproximal spaces, making it easier to remove
- Dispensed with little effort, saving the user time





Total Set Time: 2 Min.



Working Time

Futar® The comfortable classic

- Shore-A hardness of 90
- Thixotropic and easy to remove
- Sets quickly, can also be milled and easily cut with a scalpel



Total Set Time: 2 Min.



Working Time

Futar® D Hard can be so easy

- Ultra rigid Shore-D hardness of 43
- Cures with a high-luster finish, has excellent permanent dimensional stability
- Rated 4-stars in 2005 by REALITY
- Low elasticity minimizes rebound when articulating two models



Total Set Time: 1 Min.



Working Time

Futar® D Fast Twice as fast

- Rigid Shore-D hardness of 43
- Dispensed with very little effort, saving time
- Thixotropic and easy to remove



Total Set Time: 4 Min. 30sec



Working Time

Futar® D Slow Take your time

- Shore-D hardness of 43
- An exceptional combination of hardness and setting characteristics
- Ideal for myocentric bite registration, fixation of implant posts or custom impression peripheries
- Additional applications for dental laboratories



Total Set Time: 1 Min.



Working Time

Futar® Scan Perfect for CAD/CAM dentistry

- Shore-D hardness of 35
- Exceptional flexural stability
- Can be used with CAD/CAM systems
- Excellent optical properties, high recording quality and optimal image replication



Kettenbach LP offers a complete line of high-quality Futar® brand bite registrations materials that saves the practitioner time and money, without sacrificing quality. Futar® brand bite registrations will achieve proper centric relation/centric occlusion registration and minimize any adjustments.

The right materials for the right procedures!

Futar® Brand Bite Registrations

Product Properties	Futar®	Futar® D	Futar® D Fast	Futar® D Slow	Futar® Scan
Hardness	Shore-A 90	Shore-D 43	Shore-D 43	Shore-D 43	Shore-D 35
Total Working Time _{23 C / 74 F}	30 Sec.	30 Sec.	15 Sec.	1 Min. 30 Sec.	15 Sec.
Intraoral Setting Time _{35 C / 95 F}	1 Min. 30 Sec.	1 Min. 30 Sec.	45 Sec.	3 Min.	45 Sec.
Total Set Time	2 Min.	2 Min.	1 Min.	4 Min 30 Sec.	1 Min.

Available to order in normal packs

Contents	Futar®	Futar® D	Futar® D Fast	Futar® D Slow	Futar® Scan
2 x 50ml, 12 mixing tips	11912	11932	11961	11951	11971



www.kettenbach.com/www.kettenbachusa.com
 Canadian orders, please contact Cerum at 1-800-661-9567
 For U.S. orders buy direct from Kettenbach at:
 7777 Center Avenue, Suite 280 • Huntington Beach, CA 92647
 Toll free fax: 1-888-892-9820 • Local fax: 714-891-7860
 Toll free phone: **1-877-532-2123** • Local phone: 714-892-9820