

Panasil®

Impression materials designed especially
for transfer impressions



Panasil® impression materials are designed to work exceptionally well for “pick up” impressions including impression transfer impressions. Panasil® impression materials combine ideal final hardness, high tear resistance and excellent recovery from deformation. Whether you are taking an open tray transfer or closed tray transfer, you can rest assured your impressions will precisely replicate the oral situation required for precise fitting restorations such as implant restorations.

Panasil® monophas medium

This material is designed especially for transfer impressions. Thixotropic enough to build pressure in the tray, yet flowable under pressure to syringe and flow around impression copings.




Panasil® tray soft heavy

For those who would like a lower viscosity material, we offer Panasil tray soft heavy. Heavy enough to add extra pressure, yet soft enough to deform and rebound without tearing.



Panasil® initial contact

Excellent wash materials for fine detail reproduction. These materials are extremely hydrophilic and thixotropic.



Introductory offer! Purchase

2 normal packs of Panasil® tray soft heavy 50 mL cartridges
or
2 normal packs of Panasil® monophas medium 50 mL cartridges
**And you will receive one free box of our wash material
or bite registration, your choice!**
1 normal pack Panasil® initial contact 50 mL cartridges
or
1 normal pack Futar® D 50 mL cartridges

Medium Body

- 11846 Panasil® monophas medium
Normal pack cartridges 2 X 50 mL, 6 mix tips
- 28280 Panasil® monophas medium
Intro pack for Plug & Press® Dispenser
1 X 362 mL, 1 cartridge body, 10 dynamic mixers
- 11881 Panasil® monophas medium
Refill pack for Plug & Press® Dispenser
2 X 362 mL




Heavy Body

- 13541 Panasil® tray soft heavy
Normal pack cartridges 2 X 50 mL, 6 mix tips
- 28102 Intro pack for Plug & Press® Dispenser
1 X 362 mL, 1 cartridge body, 10 dynamic mixers
- 11641 Intro pack for Plug & Press® Dispenser
1 X 362 mL, 1 cartridge body, 10 dynamic mixers

Wash Materials

- 13401 Panasil® initial contact X-light
Normal pack cartridges 2 X 50 mL, 6 mix tips
- 13411 Panasil® initial contact light
Normal pack cartridges 2 X 50 mL, 6 mix tips
- 13431 Panasil® initial contact regular
Normal pack cartridges 2 X 50 mL, 6 mix tips

Accessories

- 14101 Panasil® adhesive, 10 mL 
- 17253 Dynamic mixers, 50 pieces 
- 17744 Cartridge body 
- 17234 Mixing tips, teal bag of 50 
- 17221 Intraoral tips, clear bag of 96 
- 17240 Mixing tips, yellow bag of 100 
- 17222 Intraoral tips, yellow bag of 96 



www.kettenbach.com/www.kettenbachusa.com

Canadian orders, please contact Cerum at 1-800-661-9567
For U.S. orders buy direct from Kettenbach at:
7777 Center Avenue, Suite 280 • Huntington Beach, CA 92647
Toll free fax: 1-888-892-9820 • Local fax: 714-891-7860
Toll free phone: 1-877-532-2123 • Local phone: 714-892-9820