

Kettosil®

Polyvinyl silioxane anatomical model impression material

Kettosil combines the advantages of A-Silicones and alginates to link low tear resistance with high dimensional stability. You can pour these A-Silicone impressions more than just once. Kettosil is ideal for taking impressions for:

Use as a matrix for fabrication of crown and bridge provisionals.

- Removable restorations
- Nightguards and bite splints
- Sports mouthguards
- Study Models
- Orthodontic appliances

Kettosil has many advantages over using conventional alginates:

- Multiple pours
- Pour within 30 minutes to 2 weeks
- Dimensionally stable
- Tear resistant



Kettosil® — The advantage at a glance:

- Due to its high dimensional stability, the anatomical impression can be kept for weeks and, for example, be used again if a temporary restoration is fractured or lost.
- The low tear-resistance prevents the restoration or orthodontic appliance being dislodged inadvertently, facilitates removing the impression from the mouth and reduces the risk of spontaneous extraction under the appropriate clinical conditions.
- The good thixotropic properties prevent the material inadvertently running into the pharyngeal cavity yet permit it to flow properly under pressure.
- The elastic recovery is perfectly adequate and permits the impression to be re-cast. It can also be used for impressions for straightforward removable restorations.
- The convenient foil bag system (Plug & Press) guarantees simple handling.
- The setting characteristics are pleasant for both operators and patients – the adequate total working time and brief waiting period while the impression cures in the mouth guarantee that the operator can work calmly.



Kettosil®

Medium viscosity, fast setting polyvinyl siloxane anatomical impression material.

Each pack of Kettosil includes 5 free Triple Tray impression trays for use as a matrix when fabricating provisionals and for taking quadrant anatomical impressions.

Kettosil® intro kit for Plug & Press Dispenser

28096 1 x 362mL, 10 dynamic mixers, 1 cartridge body



Kettosil® refill kit for Plug & Press Dispenser

13801 2 X 362mL



Accessories

Dynamic Mixers

17253 50 pieces



Cartridge body

17744 1 piece



www.kettenbach.com/www.kettenbachusa.com

Canadian orders, please contact Cerum at 1-800-661-9567

For U.S. orders buy direct from Kettenbach at:

7777 Center Avenue, Suite 280 • Huntington Beach, CA 92647

Toll free fax: 1-888-892-9820 • Local fax: 714-891-7860

Toll free phone: **1-877-532-2123** • Local phone: 714-892-9820