

Carriere LX™ Passive Self-Ligating Bracket System




- Elegantly Simple Locking Mechanism Opens With an Explorer and Closes Securely With the Touch of Your Finger
- One of the Lowest Profile Self-Ligating Brackets Available for Increased Patient Comfort
- Offers Quick and Easy Archwire Changes
- Low Friction, Low Force Orthodontics
- Faster Treatment Time
- Nickel-Free for Increased Biocompatibility and Acceptance
- Compound-Contour/Torque-In-Base



Ortho Organizers Carriere LX™ Process Features

**Micro-etched surface –
for increased bond strength**



- **Colour-coded for easy identification and fast placement**
- **Rigorous quality control to ensure product efficiency**










Recommended Archwire Sequence: All Archwires Pro Form™ Arch Shape

	Archwire	Package	Catalogue Numbers	
			Upper	Lower
Initial	.014 Bio-Kinetix™ with dimple	10/pk	100-852DM	100-853DM
	.016 Bio-Kinetix™ with dimple	10/pk	100-854DM	100-855DM
Interactive	.014 x .025 Bio-Kinetix™ with dimple	10/pk	100-872DM	100-873DM
	.016 x .025 Bio-Kinetix™ with dimple	10/pk	100-874DM	100-875DM
	.018 x .025 Bio-Kinetix™ with dimple	10/pk	100-866DM	100-867DM
Active	.017 x .025 CNA	5/pk	100-945	100-955
	.019 x .025 CNA	5/pk	100-948	100-958
	or	or		
	.017 x .025 Multi-Force	10/pk	100-834	100-835
	.019 x .025 Multi-Force	10/pk	100-838	100-839

Carriere LX™ Passive Self-Ligating Bracket System (continued)

McLaughlin, Bennett and Trevisi* Prescription









		Description - Colour Code	Torque	Angulation	Catalogue Numbers	
					.018 MC	.022 MC
Brackets**						
		Carriere LX Quick Pac No hooks (5x5)			703-297MQ	703-397MQ
		Carriere LX Quick Pac Hooks on cuspids (5x5)			703-298MQ	703-398MQ
		Carriere LX Quick Pac Hooks on cuspids and bicuspids (5x5)			703-299 MQ	703-399MQ
Centrals		Upper Right - Black	+17°	+4°	703-209	703-309
		Upper Left - Black	+17°	+4°	703-210	703-310
Laterals		Upper Right - Pink	+10°	+8°	703-227	703-327
		Upper Left - Pink	+10°	+8°	703-228	703-328
Lower Anteriors		Lower Right - Yellow	-6°	0°	703-237	703-337
		Lower Left - Yellow	-6°	0°	703-238	703-338
Cuspids		Upper Right - Green	-7°	+8°	703-253	703-353
		Upper Right/hook - Green	-7°	+8°	703-253HK	703-353HK
		Upper Left - Green	-7°	+8°	703-254	703-354
		Upper Left/hook - Green	-7°	+8°	703-254HK	703-354HK
		Lower Right - Blue	-6°	+3°	703-271	703-371
		Lower Right/hook - Blue	-6°	+3°	703-271HK	703-371HK
		Lower Left - Blue	-6°	+3°	703-272	703-372
		Lower Left/hook - Blue	-6°	+3°	703-272HK	703-372HK
Bicuspids		Upper Right - Purple	-7°	0°	703-259	703-359
		Upper Right/hook - Purple	-7°	0°	703-259HK	703-359HK
		Upper Left - Purple	-7°	0°	703-260	703-360
		Upper Left/hook - Purple	-7°	0°	703-260HK	703-360HK
		Lower Right - Gray	-12°	0°	703-287	703-387
		Lower Right/hook - Gray	-12°	0°	703-287HK	703-387HK
		Lower Left - Gray	-12°	0°	703-288	703-388
		Lower Left/hook - Gray	-12°	0°	703-288HK	703-388HK
		Lower Right - White	-17°	0°	703-289	703-389
		Lower Right/hook - White	-17°	0°	703-289HK	703-389HK
		Lower Left - White	-17°	0°	703-290	703-390
		Lower Left/hook - White	-17°	0°	703-290HK	703-390HK

*Ortho Organizers prescription does not imply an endorsement by the doctor(s).

**Call for availability.

Carriere LX™ Passive Self-Ligating Bracket System (continued)

Roth* Straight Arch Prescription

		Description - Colour Code	Torque	Angulation	Catalogue Numbers	
					.018 Roth*	.022 Roth*
Brackets**						
		Carriere LX Quick Pac No hooks (5X5)			703-297RQ	703-397RQ
		Carriere LX Quick Pac Hooks on cuspids (5X5)			703-298RQ	703-398RQ
		Carriere LX Quick Pac Hooks on cuspids and bicuspids (5X5)			703-299RQ	703-399RQ
Centrals		Upper Right - Black	+12°	+5°	703-207	703-307
		Upper Left - Black	+12°	+5°	703-208	703-308
Laterals		Upper Right - Pink	+8°	+9°	703-223	703-323
		Upper Left - Pink	+8°	+9°	703-224	703-324
Lower Anteriors		Lower Right - Yellow	0°	0°	703-235	703-335
		Lower Left - Yellow	0°	0°	703-236	703-336
Cuspids		Upper Right - Green	-2°	+11°	703-251	703-351
		Upper Right /hook - Green	-2°	+11°	703-251HK	703-351HK
		Upper Left - Green	-2°	+11°	703-252	703-352
		Upper Left /hook - Green	-2°	+11°	703-252HK	703-352HK
		Lower Right - Blue	-11°	+5°	703-267	703-367
		Lower Right /hook - Blue	-11°	+5°	703-267HK	703-367HK
		Lower Left - Blue	-11°	+5°	703-268	703-368
		Lower Left /hook - Blue	-11°	+5°	703-268HK	703-368HK
Bicuspids		Upper Right - Purple	-7°	0°	703-259	703-359
		Upper Right /hook - Purple	-7°	0°	703-259HK	703-359HK
		Upper Left - Purple	-7°	0°	703-260	703-360
		Upper Left /hook - Purple	-7°	0°	703-260HK	703-360HK
		Lower 4 Right - White	-17°	0°	703-289	703-389
		Lower 4 Right /hook - White	-17°	0°	703-289HK	703-389HK
		Lower 4 Left - White	-17°	0°	703-290	703-390
		Lower 4 Left /hook - White	-17°	0°	703-290HK	703-390HK
		Lower 5 Right - Red	-22°	0°	703-293	703-393
		Lower 5 Right /hook - Red	-22°	0°	703-293HK	703-393HK
		Lower 5 Left - Red	-22°	0°	703-294	703-394
		Lower 5 Left /hook - Red	-22°	0°	703-294HK	703-394HK

*Ortho Organizers prescription does not imply an endorsement by the doctor(s).

** Call for availability.