


NEW!
Exclusively from
Aurum Ceramic / Classic!


Contessa™ Zirconia

Created with Zeno® Tec for optimal anterior or posterior all-ceramic restorations



Precise internal fit ensures optimal pressure distribution.

Accurate, life-like all-ceramic results anywhere in the mouth!




Framework on model clearly illustrating 1 mm cutback for pressed ceramic.

- Traditional conservative 360° butt or chamfer margin preparations.
- Copings stopped one mm short of margin. Porcelain then pressed from margin to incisal edge for unsurpassed esthetics.
- **Use your favourite crown and bridge cement!**

Proven Strength. Superb Esthetics.

- Frameworks precision CAD/CAM milled using Yttria-stabilized Zirconium – then “hard sintered” allowing high (1300 MPa) flexural strength, reliability, fit and clinical performance.
- Larger Zirconia disc size allows milling of more extensive substructures than possible with competing systems.

For full details on exceptionally life-like, strong and durable Contessa Zirconia crowns and bridges, call your closest Aurum Ceramic/Classic laboratory.



CALL TOLL FREE

Calgary 1-800-661-1169	Edmonton 1-800-661-2745	Saskatoon 1-800-665-8815	Vancouver 1-800-663-1721	Victoria 1-800-663-6364	Kelowna 1-800-667-4146	Vernon 1-800-663-5413	Ottawa 1-800-267-7040	Toronto 1-800-268-4294
----------------------------------	-----------------------------------	------------------------------------	------------------------------------	-----------------------------------	----------------------------------	---------------------------------	---------------------------------	----------------------------------