

The Clinically Proven Answer For Snoring and Mild Obstructive Sleep Apnea

Approximately 50 million people across North America annually report some type of sleep problem including snoring and the symptoms of Obstructive Sleep Apnea (OSA). Now, with the Klearway™ appliance, available exclusively in Canada from Space Maintainers Laboratories Canada, you can offer the clinically proven solution.

Proven Through Clinical Trials

To date, more than 10,000 patients worldwide have already benefited from the use of a Klearway™ appliance. Developed by Alan Lowe, DMD, PhD, FRCD(C), FACD, Professor and Chair of the Division of Orthodontics at the University of British Columbia's Faculty of Dentistry, Klearway™ is the most thoroughly field tested and researched oral appliance for the treatment of snoring and OSA available.

Progressive titration capabilities together with excellent compliance give Klearway™ the highest reported treatment success of any oral appliance used for the treatment of snoring and OSA:

- Successfully treats snoring and 80% of mild to moderate OSA patients.
- Successfully treats snoring and 61% of severe OSA patients.
- Increases airway size in OSA patients particularly at the level of the velopharynx.
- Covert compliance monitoring shows that Klearway™ is worn on average 6.8 hours per night.
- Can be used to reduce vertical jaw posture which is significantly higher in OSA patients.



Custom Fabricated For Maximum Patient Comfort

- Created with thermoactive acrylic.
- Once warmed under hot water for easy insertion, the acrylic resin contracts as it cools, providing superior retention and fit.
- Significantly decreases soft tissue and tooth discomfort.
- Specifically designed not to encroach on tongue space.

Precisely Adjustable

- Utilizes a laser-welded (not soldered) advancement screw for superior longevity.
- Permits patient to adjust appliance according to his or her own comfort level under the guidance of the attending dentist.
- Total of 44 positions available in increments of 0.25 mm (full 11 mm range).
- Facilitates slow and gradual forward movement of the mandible avoiding patient discomfort.
- Allows patients to yawn, cough, swallow and drink water without dislodging the appliance.
- Useful for Bruxism and TMJ patients.

A Medical/Dental Team Effort

Because of the obvious life-threatening implications of sleep disorders, it is imperative that treatment commence only after the potential presence of other associated medical conditions, including sleep apnea, are ascertained. Prior to treatment with an oral appliance, the patient should be assessed by an attending physician and/or sleep specialist.

For more information on the KLEARWAY™ Appliance, call your closest Aurum/Classic - Space Maintainers laboratory TOLL FREE.



Utilizes a laser-welded (not soldered) advancement screw for superior longevity.

SPACE MAINTAINERS[®]
LABORATORIES - CANADA

© Space Maintainers Laboratories Canada (2008).
Reproduction of this work in whole or in part by any means whatsoever is strictly prohibited without the express written consent of Space Maintainers Laboratories Canada. All rights reserved.

† Some restrictions may apply due to the regulations of the College of Dental Technologists of Ontario.

Canadian Patent #2,181,360.