

# The Diagnostic Wax-up

## Show your patients what Cosmetic Dentistry can achieve!

Esthetic dentistry is, by its very nature, a very visual endeavour. Diagnostic Wax-ups, part of Aurum Ceramic/Classics' exclusive **ACCES™ (Advanced Cosmetic Communication and Esthetic System)**, create a clear blueprint for unsurpassed esthetics and functional results.

### Powerful Patient Communication Tool

- Natural-looking, three-dimensional representation of final case.
- Illustrates current situation and how more esthetic smile can be achieved.
- Shape and contour of the teeth easily discussed with patient.
- Dramatically increases treatment acceptance rates.

### Visual Aid for Adequate Tooth Preparation

- Obtain a visual understanding of tooth reduction requirements, reducing preparation appointment time in the operator.
- Matrices serve as a valuable check to ensure adequate preparation.
- Pre-plan location and contour of any gingival tissue recontouring.

### Template for Final Restorations

- Assists in creating a smile design that assures predictability in the mouth.
- Allows careful evaluation of function, anterior guidance and cuspid rise.
- Provides a template from which both the preparation and the temporary index can be fabricated.

### Allows Quick Fabrication of Outstanding Temporaries

- Aurum Ceramic/Classic will provide a Temporary Stent from the Diagnostic Wax-up.
- Temporaries that mimic the final restorations can be created within ten minutes.
- Allows the patient a few days to evaluate the shape and function of their "new teeth".

### Build Your Practice With Aurum Ceramic/Classic

For more information on ACCES™ and Diagnostic Wax-ups, call your closest Aurum Ceramic/Classic laboratory TOLL FREE.



### Case Study (Dentistry courtesy of Dr. Orin Harker)



*Patient smile before preparation.*



*Diagnostic Wax-up depicting natural-looking, three-dimensional representation of the final case.*



*One-piece Temporary made from Diagnostic Wax-up and Temporary Stent. No bonding required. Temporization completed in less than 10 minutes.*



*Final IPS Empress restorations duplicating length and shape of temporaries. Illustrates the outstanding results possible with ACCES, Diagnostic Wax-ups and pre-treatment planning.*