



The New Era In Splint Therapy

CLEAR-FLEX

Orthodontic splints have always allowed the practitioner a trial phase of therapy in derangement cases, TMJ and other craniomandibular disorders, without making irreversible changes in condylar position. Unfortunately, splints were typically constructed from a hard acrylic material — a material that was subject to breakage, difficult to work with and often uncomfortable to wear. Now, with CLEAR-FLEX splints from Space Maintainers, a new era in splint therapy has arrived!

Unique Clear Thermo-Elastic Material With Memory

- resists breakage
- allows splint to maintain original shape
- becomes flexible when placed under hot tap water
- as it cools, splint returns to its original shape and conforms to teeth
- greatly improved patient comfort, acceptance and compliance

Wide Range of Application

- reproduce mandibular rest position as a centric occlusion position
- provide therapeutic occlusion in harmony with the neuromuscular patterns
- provide disclusion guidance which frees masticatory system from unwanted posterior interferences
- provide a masticatory machine which can apply maximum force with minimum self destruction and pain

Very Accurate Fit

- chairtime drastically reduced
- occlusal adjustments, if any, done in usual manner

Clasping For Retention Almost Eliminated

- usually only in cases where there is no height of contour
- ball clasps can be added easily on linguals and buccals of the bicuspid and molars

Simple Prescription Requirements

- upper and lower casts and bite registration (or indication of the vertical) all that is required

Experience The New Era In Splint Therapy!

For more information on CLEAR-FLEX SPLINTS and how Space Maintainers can help you increase your practice, call TOLL FREE.

SPACE MAINTAINERS™
LABORATORIES - PACIFIC NORTHWEST